

# Quorum

## Qdry



## Automated Critical Point Dryer

### Features and Benefits

- ◆ Automated control offering repeatable outcomes with minimal user intervention
- ◆ Built-in adiabatic cooling for fast cool down rates allowing for more throughput
- ◆ Easy-to-use and flexible software allowing for quick set-up of profiles
- ◆ Status display to indicate process time remaining
- ◆ Process driven, with individual user profiles for quick set-up time
- ◆ Supplied with pre-set profiles covering different sample types
- ◆ Screw top chamber for easy sample loading
- ◆ Large viewing window ensuring observation of process flow
- ◆ Small footprint ideal for use in fume hoods



**SPECIALISTS IN EM SAMPLE PREPARATION**

**E: SALES@QUORUMTECH.COM • T: +44 (0) 1323 810981 • W: QUORUMTECH.COM**

# Qdry



## Specifications

### Instrument Dimensions

◆ W 450 mm x H 325 mm x D 505 mm

### Weight

◆ 35 kg

### Specimen holders

◆ AL800019-1 specimen holder as standard

### Safety

◆ Built-in safety devices in both hardware and software, including relief valves, level sensors, pressure switches and programming restraints

### Chamber diameter

◆ 60 mm Ø

### Internal chamber height

◆ 59 mm, usable chamber height = 34.5 mm, usable chamber volume = 97.5 ml




### Critical point dryer liquid carbon dioxide (CO<sub>2</sub>) requirements

◆ All models require a cylinder of liquid CO<sub>2</sub> fitted with a siphon tube (normally indicated by a vertical white stripe on the cylinder). If there is any doubt regarding the presence of a siphon tube, advice should be sought from the supplier.

**Please contact our sales team at [sales@quorumtech.com](mailto:sales@quorumtech.com) for full pricing.**

**SPECIALISTS IN EM SAMPLE PREPARATION**

**E: SALES@QUORUMTECH.COM • T: +44 (0) 1323 810981 • W: QUORUMTECH.COM**





## HEADQUARTERS

### Quorum Technologies

Judges House  
Lewes Road  
Laughton  
East Sussex  
BN8 6BN  
United Kingdom

**T:** +44 1323 810981


**E:** [sales@quorumtech.com](mailto:sales@quorumtech.com)

**W:** [www.quorumtech.com](http://www.quorumtech.com)

**Please contact our sales team at [sales@quorumtech.com](mailto:sales@quorumtech.com) for full pricing.**

**Products are for research use only.**

While every effort is made to ensure accuracy, no responsibility can be accepted for errors or omissions. Data may change, as well as legislation and you are strongly advised to obtain copies of the most recently issued regulations, standards and guidelines. This document does not form the basis of a contract. Quorum has a policy of ongoing product improvement and reserves the right to change specifications without notice © 2023 Quorum Technologies Limited. All rights reserved.



**SPECIALISTS IN EM SAMPLE PREPARATION**

**E: SALES@QUORUMTECH.COM • T: +44 (0) 1323 810981 • W: QUORUMTECH.COM**